

Reborn as Innovation

# Raphael 606-PR

R606-P

In 2015 January, 01M signed an Asian license for manufacturing & distribution from American Leander Tables. Since 2016, 01M developed first 606 Auto Flexion table and exported to China, Europe, Taiwan and many other countries. In 2018 November, 01M hosted the Leander Technique Seminar from Dr. Dennis Anthony, CEO of Leander Tables. In 2019 February, the latest new 4 class Raphael 606 tables were launched.



Color



Specifications

Width	640 mm
Length	1,860 mm ~ 2,260 mm
Height	585 mm ~ 735 mm
Weight	142 kg

World Class Spine Care Table R606

Advanced auto flexion table made by clinical studies

US FDA Listed

It is a table applied with auto-flexion, auto-drop, flexion memory, advanced scoliosis treatment kit and IT technology.

Raphael 606  
Story of Development Design

Raphael 606 - PR is the first Korean Leander table which was developed in the United States through a licensing contract. It was improved with clinical studies from doctors and technological innovation of YOUNG IL M.



Registered Patent of Design

Feature



Center-Forward Auto Drop

YOUNG ILM's patented technology  
Center-Forward vertical drops applied



Abdominal cushion adjustable

Break-Away abdominal cushion for pregnant, large abdomen patients



Pattern and memory function

Easy use with pattern and memory functions with YOUNG ILM's technology



Manual Flexion Handle

Left and right lateral bending is possible through the flexion rod



Integrated Electronic Control Panel


Controllable the drop and flexion at the same time with integrated electronic control panel

Option



Latest Scoliosis Belt

The set of scoliosis belt is used to treat scoliosis patients




Auto-Flexion, Flexion Belt


Increase the effect of flexion by fixing the waist



Auto Drop



Auto Flexion



Lateral Bending



Break Away